



X⁶ Deluxe Chair Bed



X⁶ Deluxe Chair Bed

The Care Quip X⁶ Deluxe Chair Bed has been designed for the treatment and prevention of pressure problems related to long-term seating.

Using a contoured dual layer air cushion overlaid with pressure relieving high density foam, the X⁶ Chair Bed evenly distributes the occupant's weight, ensuring that there are no main pressure points that can restrict circulation.

The Care Quip X⁶ Deluxe Chair Bed is a must have in facility prevention and treatment of pressure sores. This cost effective, mobile chair bed improves facility care and equips carers with quality pressure care options.



Seating Surface

Four way Stretch Vapour permeable fabric Seating surface increases patient comfort and is designed to improve patient pressure care.

Vapour permeable - water resistant yet breathable and improves skin microclimate

Belgium Made fabric meets European fire standards for being Halogen free fire retardancy

All covers are removable for infection control management

Pressure Relieving Foam Back

Shaped foam and air chambers promotes exceptional pressure redistribution

Air flow chambers helps reduce skin moisture when seated for extended periods of time

High Density foam outer supports improves back posture and redistributes pressure

Foam and air chambers can be removed and customised to the patient's pressure requirements

Push Handle

Full Hand tilt and Recline levers allows seamless Patient comfort

Vertical and Horizontal Push Handle positions increases movability and chair control



Foot Plates

Removable comfort padded foot plates to aid patient comfort and heel pressure remediation

Safe Handling

Three individual locking castors, along with one directional castor assists with manoeuvring, and the ergonomic handles assist with leverage

Swing Away Backrest Laterals and Armrests

Increases patient access, and decreases carers demands when positioning the patient

Allows for easy patient side transfers

Air Seat Section

Contoured Adjustable air cushion increases patient emersion and pressure care

Air Cushion is overlaid with a high-density pressure relieving foam that increases patient stability and comfort

Specifications

Seat Height	650mm
Seat Depth	560mm
Seat Width	500mm
Seat to Top of Back	730mm
Overall Height	1290mm
Overall Width	700mm

Overall Depth	1130mm
Recliner Length	1790mm
Tilt Angle	25°
Leg Raise Angle	90°
Recliner Angle	50°
S.W.L	180KG