

# New Product

## DXR4000 Steiss Quartet 4 in1 System



Grade 4\*

Capacity  
220kg

Cycle  
20min

Ratio  
1:4

Hour  
Counter

### Overview

Product Code: MA0600

The Steiss™ DXR4000 Quartet is a high specification Alternating Air-Pressure Full Mattress Replacement System and is designed to prevent and treat pressure ulcers up to grade IV.

The system uses the latest technology in alternating mattress therapy to provide accurate pressure settings for individual patients' needs. The 20 cells of the mattress are arranged in a transverse manner to provide extra protection and comfort for the patient. Through the use of alternating therapy; Steiss™ DXR4000 Quartet system works on the principle of actively encouraging tissue blood flow by gently increasing and decreasing pressure within the support surface over the 20 minute cycles.

This results in the reduction of surface pressures helping to prevent tissue breakdown and encourage healing. For added durability, and longer working life, your Steiss™ DXR4000 Quartet mattress has been constructed using high strength components and premium grade materials.

The mattress features a heavy duty base layer and a 100% waterproof, multi-stretch vapor permeable PU cover to reduce risks of friction and shear.

This 'Grade' is a guide only

Each Individual must be assessed by a tissue viability specialist, the correct product chosen and set to operate according to the patient's needs.

Regular repositioning schedules need to be in place as this product does not replace good nursing care.

Continue to next page...

# New Product

## DXR4000 Steiss Quartet 4 in1 System



### Features

- Digital Alternating Air-Pressure Full Mattress Replacement System, suitable for home and care units
- Grade IV risk / management\*
- 1:4 alternation therapy in 20 minute cycles
- 20 supportive, full depth cells for superior pressure care management
- 3 static head cells for greater patient comfort
- 7 specialist narrow heel cells for targeted pressure care
- Separate pocket, for optional foam insert, underneath air cells
- Static and alternating modes
- Figure-of-8 air circulation
- Static mode with 30 minute auto-return for safer patient transfers and nursing tasks
- 6 clear and intuitive digital adjustment comfort settings
- Auto-firm function providing faster static inflation for aiding nursing tasks
- Multi stretch, highly vapour-permeable PU top cover providing low shear, friction and moisture protection
- External CPR valve design for fast emergency deflation
- Premium quality all round zip with large waterfall cover for improved infection control
- Strong base-mount straps holding the mattress when insitu and when in the rolled-form during storage
- Anti-kink tubing for added safety and reliability
- Simple filter cleaning and exchange
- Low Pressure, Power Failure & Service/Maintenance Due audio-visual alarms
- Auto lock-out feature
- Transport Security Device
- 3 year guarantee
- CE marking
- Mattress size 200 x 90 x 24cm (79" 35.5 x 9.5")
- Minimum patient weight 40kg/6.3 stone
- Maximum patient weight 220kg/34.64 stone

---

Care Quip Pty Ltd  
5 Millennium Court Silverwater NSW 2128 Australia  
P (02) 9714 1111 F (02) 9714 1133 E sales@carequip.com.au  
www.carequip.com.au