





Designed for Comfort. Built for Care.

The Movi® Contour X4 is a therapeutic lift-and-recline chair designed to support individuals with reduced mobility, postural instability, fatigue, and functional transfer limitations in domestic and supported living environments.

Incorporating tilt-in-space positioning, the Contour X4 enables effective pressure redistribution, pelvic stability, and postural alignment, reducing shear forces and supporting prolonged sitting tolerance. This makes the chair suitable for participants requiring postural support, fatigue management, and pressure care intervention as part of an ongoing seating prescription.

The chair's smooth lift-assist and recline mechanisms facilitate safe sit-to-stand transfers, promoting user independence while reducing manual handling risk for carers and support workers. Controlled movement through all positioning ranges ensures stability and confidence for users with balance or strength limitations.

A modular comfort and support system allows the Contour X4 to be configured to individual clinical needs, accommodating changes in functional ability over time. High-quality, contoured cushioning supports comfort, trunk positioning, and long-term durability in high-use care settings.

Suitable for use across NDIS-funded supports, home care packages, and residential aged care, the Movi® Contour X4 provides a clinically appropriate, prescription-ready seating solution that supports independence, safety, and quality of life.

Features:

- » Lift, Recline & Tilt-in-Space Functionality
- » Independently controlled backrest and footrest for personalised comfort
- » Tilt-in-space positioning evenly distributes body weight
- » Helps reduce spinal pressure, improve circulation, and relieve lower-back discomfort
- » Adjustable two position backrest improves seat depth when needed
- » Adjustable backrest cushions providing enhanced client pressure care and support



Code	Size	Color	Motor	Tilt in Space Function	Removable Seat Foam	Seat Height	Seat Depth (Position 1)	Seat Depth (Position 2)	Seat Width	Seat to Top of Back	Max User Weight
ED5340	Small	Black	Dual Motor	Yes	Yes	450mm	450mm	500mm	460mm	700mm	180kg
ED5350	Medium	Black	Dual Motor	Yes	Yes	470mm	450mm	500mm	500mm	700mm	180kg
ED5360	Large	Black	Dual Motor	Yes	Yes	500mm	450mm	500mm	540mm	700mm	200kg
ED5310	Small	Beige	Dual Motor	Yes	Yes	450mm	450mm	500mm	460mm	700mm	180kg
ED5370	Medium	Beige	Dual Motor	Yes	Yes	470mm	450mm	500mm	500mm	700mm	180kg
ED5330	Large	Beige	Dual Motor	Yes	Yes	500mm	450mm	500mm	540mm	700mm	200kg



Electric Hand Control

With a built-in USB charging port, featuring simple, clearly marked buttons for smooth and safe movement, plus a quick-release handset to prevent unauthorized use.



Modular Backrest System

Three detachable back cushions (top, middle, bottom) secured by hook and loop fasteners, allowing easy removal, cleaning, replacement, or comfort adjustment.



Adjustable Seat Depth

Suit different user heights, ensuring optimal posture, comfort, and support.



Designed for Safety & Stability

Smooth, controlled movement for safe lifting and reclining, hazard-aware construction to minimize risks, and a stable, supportive frame for peace of mind.



Three size options for better support.

Available in small, medium, and large sizes with weight capacities up to 200 kg, ensuring the right fit for every user.



Durable Upholstery Options

Available in easy-care fabric or leather, designed for everyday use with low-maintenance cleaning and long-lasting appearance.



Product Warranty

5 year warranty across the complete chair, providing peace of mind and assurance to every client

Experience comfort without compromise.

Discover how the Movi® Contour X4 Recliner can transform everyday seating into a safer, more comfortable experience.

1300 930 930

Contact us today to learn more, request a quote, or arrange a consultation.

