





Designed for Comfort. Built for Care.

The Movi® Recline X3 is a premium lift-and-recline chair designed to deliver exceptional comfort, safety, and support for everyday living. Engineered with advanced tilt-in-space technology and smooth electric controls, it helps users sit, recline, and stand with ease—promoting independence, comfort, and well-being at home.

Whether relaxing, reading, or transitioning safely from seated to standing positions, the Movi® Recline X3 offers reliable support without compromising on style or durability.

Features:

- » Independent control of the backrest and footrest
- » Tilt-in-space positioning evenly distributes body weight
- » Reduces spinal pressure and lower-back discomfort
- » Improves circulation and long-term seating comfort
- » Sit to stand Functionality



Code	Size	Color	Motors	Tilt in Space Function	Seat Height	Seat Depth	Seat Width	Seat to Top of Back	Max User Weight
ED5270	Small	Ash	Dual Motor	Yes	470mm	450mm	460mm	600mm	180kg
ED5200	Medium	Ash	Dual Motor	Yes	500mm	500mm	500mm	650mm	180kg
ED5290	Large	Ash	Dual Motor	Yes	530mm	520mm	540mm	700mm	200kg
ED5210	Small	Taupe	Dual Motor	Yes	470mm	450mm	460mm	600mm	180kg
ED5240	Medium	Taupe	Dual Motor	Yes	500mm	500mm	500mm	650mm	180kg
ED5250	Large	Taupe	Dual Motor	Yes	530mm	520mm	540mm	700mm	200kg



Electric Hand Control

With a built-in USB charging port, featuring simple, clearly marked buttons for smooth and safe movement, plus a quick-release handset to prevent unauthorized use.



Removable Seat Cushion with Replaceable Core

A high-density foam core delivers superior support, with Velcro attachment for easy removal. If client needs change, the seat cushion can be upgraded to an enhanced pressure care option.



Built for Safety & Stability

A hazard-aware design that minimizes trapping risks and a stable, supportive frame for safe lifting and daily use.



Three size options for better support.

Available in small, medium, and large sizes with weight capacities up to 200 kg, ensuring the right fit for every user.



Durable Upholstery Options

Easy-care fabric or leather, designed for everyday use with low-maintenance cleaning and care.



Product Warranty

5 year warranty across the complete chair, providing peace of mind and assurance to every client

Experience comfort without compromise.

Discover how the Movi® Recline X3 Recliner can transform everyday seating into a safer, more comfortable experience.

1300 930 930

Contact us today to learn more, request a quote, or arrange a consultation.

